NOVOMEDICS  
MODERNE MEDIZINTECHNIK

# Produktbeschreibungen

NOVOMEDICS GMBH STAND 30.11.2020

# Inhalt

Beckenkammkeil.....	S. 3
BioAdapt.....	S. 5
BioReady.....	S. 7
BioSet.....	S. 9
Fascia Lata.....	S. 11
Fortiva.....	S. 13
Pericardium.....	S. 15
Rippenknorpel.....	S. 17
Spongiosa.....	S. 19
Tutobone.....	S. 23
Tutomesh.....	S. 26
Tutopatch.....	S. 30

## Tutoplast® Beckenkamm Keile

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## Studien - Referenzen

Lee, Yu-Po R., Mark C. Jo, Mario Luna, Bobby Chien, Jay Lieberman, and Jeffrey Wang. "The Efficacy of Different Commercially Available Demineralized Bone Matrix Substances in an Athymic Rat Model." *Journal of Spinal Disorders & Techniques* 18.5 (2005): 439-44.

Moore, Steven T, Jordan M Katz, Rasa M Zhukauskas, RYanne M Hernandez, Christine S Lewis, Peter Robert Supronowicz, Elise Gill, Susan M Grover, Nathaniel S Long, and Ronald R Cobb. "Osteoconductivity and Osteoinductivity of Puros® DBM Putty." *Journal of Biomaterials Applications* 26.2 (2011): 151-71.

Wang, Jeffrey, C. Alanay, A. Mark, Davies Kanim, Linda Campbell, E. Dawson, and Pat Lieberman. "A Comparison of Commercially Available Demineralized Bone Matrix for Spinal Fusion." *European Spine Journal* 16.8 (2007): 1233-240.

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## Studien – Referenzen

Brescia et al., 2016. Post-incisional ventral hernia repair in patients undergoing chemotherapy: improving outcomes with biological mesh. World journal of surgical oncology, 14(1), p.257.

Ferzoco, S. et al., 2015. Mesh & Prosthesis. Hernia, 19(Supplement 1), pp.S147–S156.

Maxwell, Daniel W., Alexandra M. Hart, Orion P. Keifer, Sameer H. Halani, and Albert Losken. "A Comparison of Acellular Dermal Matrices in Abdominal Wall Reconstruction." Annals of Plastic Surgery 82.4 (2019): 435-40. Web.

Paepke, Stefan, Marion Kiechle, Christine Ankel, Joy Weyrich, Evelyn Klein, Anna Schneider, Navid Johannigmann-Malek, Anna-Sophia Dietrich, Ralf Ohlinger, and Marc Thill. "Surgical Studies of Reconstructive Breast Surgery – An Overview of the Topics at the 2019 Annual Meeting of the Working Group for Reconstructive Surgery in Oncology-Gynecology." 79.6 (2019): 584-90. Web.

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
## Studien - Referenzen

- ANAGNOSTAKOS K. UND KOCH K.  
Pharmacokinetic properties and systemic safety of vancomycin-impregnated cancellous bone grafts in the treatment of spondylodiscitis BioMed Research International; 2013  
Eight patients with a destructive spondylodiscitis (3 thoracic, 1 thoracolumbar, and 5 lumbar) have been treated with vancomycin-loaded Tutoplast cancellous Bone. Patients with large defects received a vertebral body replacement system and Tutoplast cancellous chips within. In patients with small defects an antibiotic-loaded Tutoplast femoral head was inserted in a press-fit technique instead. The results indicate that vancomycin can be locally released in high concentrations without any risk of a systemic toxicity. The infection eradication rate was 100% at a mean follow-up of 16.5 months.
- RAJAN G.P., FORNARO J., TRENTZ O. AND ZELLWEGE R.  
Cancellous allograft versus autologous bone grafting for repair of comminuted distal radius fractures: a prospective, randomized trial The Journal of Trauma 60(6): 1322-1329; 2006  
90 patients with distal radius fractures were randomly assigned to receive either autograft (n=46) or Tutoplast® Cancellous Bone (n=44). Radiological results were comparable and within normal range. Surgical time was 28 minutes longer on average in the autograft group and after one year 17% of these patients still had loss of sensitivity on the thigh or pain at the harvesting site.
- BÄRLOCHER C.B., M.D., BARTH A., KRAUSS J., BINGGELI R. AND SEILER R.W.  
Comparative evaluation of microdiscectomy only, autograft fusion, polymethylmethacrylate interposition and threaded titanium cage fusion for treatment of single-level cervical disc disease: a prospective randomised study in 125 patients.

Neurosurgical Focus, 12(1); 2002

In 36 of the 125 patients the cage was filled with Tutoplast® Cancellous Bone. This group had the highest fusion rate.

- BENLI I.T., AKALIN S., KIS M., CITAK M., AYDIN E. AND DUMAN E.  
Frontal and sagittal balance analysis of late onset idiopathic scoliosis treated with third generation instrumentation. Kobe Journal of Medical Science, Vol 47; 2001  
Report on 217 patients who had posterior only grafting with autogenous or allogenic cancellous bone in addition to internal fixation. 33 patients had Tutoplast® Cancellous Bone. The overall fusion rate was 94.5%.
- OZLU S., AKALIN S., BENLI T., KIS M., CELIK N. AND BAZ A.B.  
Radionuclid and histopathological evaluation of combined allograft and autogenous applications in anterior and posterior spinal fusion  
White Paper: Ankara Social Security Hospital; 2001, abstract: The Journal of Turkish Spinal Surgery, 8(1); 1997  
In 42 patients a combination of autogenous and Tutoplast® Cancellous Bone was used for spinal fusion. After 12 months x-rays and scintigraphy showed solid fusion in all patients. A comparative histology from 11 of these patients with 5 with autograft only showed no differences.
- GAMBINI A. MASTANTUONO M, DI GIORGIO L., VALEO M., BASSETTI E., PERSIANI P., COSTANTINO E. AND VILLANI C.  
Rehabilitation of allograft with bone dehydrated with solvents in reconstruction after removal of bone tumors: MRI evaluation Chir. Organi Mov. LXXXIV; 1999  
In 15 patients the bony defect after resection of benign or moderately malignant tumors was filled with Tutoplast® Cancellous Bone.

In 13 patients the MRI showed complete remodelling after 6 to 13 months while in 2 remodelling was only partial after 6 and 8 months.

- KIM Y.Y., SHIM J.I., SUNG Y.B. AND KIM S.K.  
Revision arthroplasty using allogenic cancellous chip with HA for pronounced acetabular bone deficiencies Bioceramics, Vol. 9; 1996  
In 100 patients with large bone defects at the acetabulum Tutoplast® Cancellous Bone augmented with hydroxyapatite granules and soaked with autogenous bone marrow was used. All but 3 showed good bony integration and remodelling after one year.

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## Studien – Referenzen

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#### Abdominal

- GURRADO A ET AL.  
Impact of pericardium bovine patch (Tutomesh®) on incisional hernia treatment in contaminated or potentially contaminated fields: retrospective comparative study. Hernia, 1; 2014

152 patients (mean age 65 ± 14 years) underwent incisional hernia repair (2002-2012) in complicated settings. 76 patients were treated with primary closure technique; 38 with reinforcement by synthetic and 38 with Tutomesh. Significant decreases in Tutomesh group were observed for wound infection, recurrence and seroma. The authors conclude that Tutomesh seems to be safe and effective to successfully repair ventral hernia in potentially contaminated operative fields, especially in association with bowel resection.

- KOZIOL D AND KAHLE M  
Biological mesh in conventional surgery of inguinal hernia: results after 6 years of use Hernia, 13 (Suppl 1); 2009  
Report on 176 patients of average age 66.3 years who received Tutomesh for inguinal hernia repair.

The collagen is typically remodeled within 6 to 8 weeks and replaced by host connective tissue. This results in a stabilization of the transversalis level and reinforcement of the deep inguinal ring.

- URBACH V ET AL.  
Bovines Perikard ein vielseitiges Material der Zukunft (Bovine Pericard, a versatile material of the future) Abstract. 126th Congress of the German Society for Surgery, Munich, Germany, 2009  
200 patients were treated between 2003 and 2008 with more than 300 pericardial membranes. Under correct application and indication Tutomesh offers a safe connective tissue substitute without the risk of unpredictable late complications and increased risk in case of reoperation.
- CAVALLARO A ET AL.  
Use of biological meshes for abdominal wall reconstruction in highly contaminated fields World J Gastroenterol, 16(15); 2010 Case Report on two patients who had a severe previous history of synthetic mesh failure following hemicolectomy. After removal of the contracted mesh, the posterior rectus sheath was repositioned to the midline, and a sheet of Tutomesh® was implanted in the retrorectus space and secured to the anterior abdominal wall in a tension-free manner. The patients had no significant postoperative complications. At 4- respectively 5-year follow-up visit, there was no evidence of recurrence. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC2856837/pdf/WJG-15-1928.pdf>
- CAVALLI M ET AL.  
Pubic Inguinal Pain Syndrome: the So-called Sports Hernia. Surg Technol Int, 24; 2014

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### Tutomesh®

<b>68443</b>	120mm x 160mm
<b>68442</b>	80mm x 120mm
<b>68444</b>	120mm x 200mm
<b>68441</b>	60mm x 80mm
<b>68440</b>	40mm x 50mm
<b>68542</b>	130mm x 220mm
<b>68541</b>	120mm x 200mm
<b>68540</b>	110mm x 180mm
<b>68543</b>	100mm x 160mm

Report on surgical treatment of chronic groin pain when conservative treatment was not effective. The surgical procedure includes the release of all three nerves of the region, the correction of the imbalance in strength with the partial tenotomy of the rectus and adductor longus muscles, and the repair of the partial weakness of the posterior wall with Tutomesh. This treatment reported excellent results with complete relief of symptoms after resumption of physical activity in all cases.

- D'AMBRA L ET AL.  
Use of bovine pericardium graft for abdominal wall reconstruction in contaminated fields  
World J Gastrointest Surg, 4(7); 2012 48 patients underwent between 2006 and 2010 abdominal wall reconstruction using Tutomesh. 34 of them had a contaminated wound and 14 hernia relapse on infected synthetic mesh. According to the authors recurrence rate was similar to that of prosthetic mesh repair but the application of Tutomesh is moreover a safe and feasible option for managing complicated abdominal wall defects where prosthetic mesh is unsuitable.  
<http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3420984/pdf/WJGS-4-171.pdf>
- FRANCONI G ET AL.  
Giant abdominal wall defects  
Crovella F. et al. (Ed.) Incisional hernia, 2008  
Giant abdominal wall defect is a complex surgical problem that must be treated in two or more steps. The standardized use of synthetic prostheses is never justified. The authors think the use of biological material can be considered the way of the future for this pathology.
- KOZIOL D ET AL.  
Biologische Netze. Ein Fortschritt in der Leistenhernienchirurgie?  
Abstract. Association of Bavarian Surgeons, 81th Annual Meeting, 2004

Tutomesh 8x6cm was used in 81 patients with direct inguinal hernias since September 2002. 58 could be followed until February 2004. No recurrences, lasting pain or foreign body feeling was seen.

- LINDEMANN M ET AL.  
Der operative Narbenhernienverschluss. Bioimplantat aus bovinem Perikard  
MedReport, 28(23); 2004  
38 patients with an incisional and 2 with a parastomal hernia were treated with Tutomesh. 1 large seroma developed that could be successfully treated. In 4 patients local pain occurred during distinct movements at the fixation suture site. No other complication occurred during the follow-up period.
- URBACH V ET AL.  
Bovines Perikard Anwendungstechniken bei der operativen Versorgung von Hernien (Bovine pericardium techniques for use in the surgical treatment of hernias)  
82/BIOmaterialien, 8(S1); 2007  
Report on the use of 300 Tutomeshes in 206 patients for the repair of 104 inguinal hernias, 82 incisional hernias and 20 other hernias over a period of 4 years. With proper size and fixation technique the recurrence rate was low. (For inguinal hernia: the adequate size of 10x15 cm and extensive fixation esp. with fibrin glue; for incisional hernia: strictly tension-free sublay technique, unsuitable for inlay technique). No complications known from synthetic meshes occurred. Biopsies showed complete remodeling after 2-3 months (in German).
- VAN TUIL C ET AL.  
Experience with management of anterior abdominal wall defects using bovine pericard  
Hernia, 10; 2006  
Report on 24 patients with congenital abdominal wall defects treated with Tutomesh and Tutopatch over a period of 5 years.

The authors conclude that Tutopatch and Tutomesh are an optimal biomaterial for the closure of anterior abdominal wall defects in gastrochisis and omphalocele.

#### Extremitäten

- NATHAN SS ET AL.  
Collagen membranes for host-implant integration: a pilot clinical study.  
J Orthop Surg, 19(2); 2011  
Retrospective study on 14 patients who underwent limb salvage surgery for musculoskeletal oncological diseases. Musculoskeletal defects were covered with Tutomesh. No patient had any infection or any foreign body reaction. The Tutomesh membrane facilitated host-implant integration and provided a feasible anatomic reconstruction for ligaments in the shoulder, knee, and hip.  
<http://www.josonline.org/pdf/v19i2p157.pdf>
- Brustrekonstruktion  
GUBITOSI A ET AL.  
Acellular bovine pericardium dermal matrix in immediate breast reconstruction after Skin Sparing Mastectomy.  
Int J Surg, 12 (Suppl 1); 2014  
A prospective study, with 24 patients who underwent 28 immediate breast reconstructions with Tutomesh (wrongly named as acellular dermal matrix ADM). The authors conclude that the use of Tutomesh for immediate breast is safe and technically useful. Complications rate is not high, except for seroma formation that can be reduced by the contemporary use of fibrin sealant.
- Tierstudien  
GUERZON ER ET AL.  
Collagen membranes for host-graft integration: an animal study.  
J Orthop Surg, 19(2); 2011

18 Achilles tendons in 9 rats were used. A 15x5 mm membrane (Tutomesch) was laid circumferentially around the tendon and secured with prolene 3-0 sutures. The collagen membrane achieved favourable host-graft integration as early as week 4. There was no host-graft reaction or wound complication (infection, abscess or seroma).

<http://www.josonline.org/pdf/v19i2p151.pdf>

▪ KAPAN S ET AL.

Comparison of PTFE, pericardium bovine and fascia lata for repair of incisional hernia in rat model, experimental study  
*Hernia*, 7(1); 2003

A 1x2cm subcutaneous full thickness abdominal wall defect was created in rats and closed with patches of Tutopatch, Tutoplast Fascia lata and Goretex Dual Mesh. After one and two weeks bursting pressure, tensile strength of the graft-tissue connection, histological graft maturation and adhesion formation were evaluated. All values were comparable for Tutopatch and GoreTex while bursting strength was lower with Fascia lata.

▪ TESTINI M ET AL.

Bovine pericardium patch wrapping intestinal anastomosis improves healing process and prevents leakage in a pig model.

*PLOS ONE*, 9(1); 2014

Use of pericardium bovine patch (Tutomesch) as reinforcement of intestinal anastomosis is safe and effective, significantly improving the healing process.

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
\*Special Products | available only on request

## Tutogen, Tutopatch - Artikelgrößen

## Studien – Referenzen

### Neurosurgery

- ARTS MP ET AL.  
Surgical treatment of idiopathic transdural spinal cord herniation: a new technique to untether the spinal cord  
Acta Neurochir (Wien), 148; 2006  
Two case reports of the successful use of Tutopatch for the closure of anterolateral defects of the spinal dura. Tutopatch was wrapped around the spinal cord and fixed with fibrin glue.
- BAUMANN I ET AL.  
Ossifying fibroma of the ethmoid involving the orbit and the skull base  
Otolaryngol Head Neck Surg, 133(1); 2005  
Case report. Tutopatch was used to reconstruct the skull base after tumor resection.
- FILIPPI R ET AL.  
Bovine pericardium for duraplasty: clinical results in 32 patients  
Neurosurg Rev, 24; 2001

Tutopatch was used for duraplasty in 32 patients. No complications attributable to the graft occurred. The authors conclude: solvent preserved bovine pericardium should be recommended as a safe, suitable, and cost-effective material for duraplasty.

- GÜRER B ET AL.  
Use of the bovine pericardial patch and fibrin sealant in meningomyelocele closure  
*Acta Neurochir (Wien)*, 156(7); 2014  
Eight meningomyelocele patients were treated with Tutopatch and fibrin sealant technique at the fascial level-between the dural sac and the skin. Stable coverage was achieved and none of the possible complications was observed postoperatively.
- MUSSACK T ET AL.  
Bilateral decompressive craniectomy due to intracranial hypertension during acute posttraumatic liver dysfunction  
*J Trauma*, 58(5); 2005  
Case report on a successful bilateral dural enlargement with Tutopatch for brain decompression.
- RINCON I ET AL.  
Reconstructive possibilities in surgery of the frontal skull base and cavity. A case report.  
*Rev Esp Cir Oral y Maxilofac*, 30(1); 2008  
Case report on the reconstruction of the frontal skull base after tumor resection. Tutopatch was used to reconstruct the dural defect.
- ROSER F AND TATAGIBA MS  
Posttraumatic Syringomyelia  
In: *Neurosurgical forum, letters to the editor; J. Neurosurg Spine*, 6(2); 2007  
A duraplasty with implantation of membrane material (e.g. Tutopatch) should be mandatory to avoid new arachnoid scarring within the 1st postoperative month and to allow restoration of an adequate cavity for cerebrospinal fluid flow.

## Ophthalmology

- ALIO, JL, ET AL.  
Bovine pericardium membrane (tutopatch) combined with solid platelet-rich plasma for the management of perforated corneal ulcers.  
*Cornea*, 32(5); 2013  
A total of 6 patients were treated with Tutopatch and eye platelet-rich plasma (E-PRP) placed on the corneal surface, underneath the patch. In all cases, the corneal perforation was sealed. No evidence of infection or inflammation was detected. Digital tonometry confirmed acceptable ocular pressure in all cases. No patients reported pain, discomfort, or any subjective symptoms, and no complications were observed.
- GUPTA M ET AL.  
Bovine pericardium (Tutopatch) wrap for hydroxyapatite implants  
*Eye*, 21(4); 2007  
Prospective study with 19 patients undergoing enucleation for uveal melanoma. The hydroxyapatite implants were wrapped with Tutopatch. Follow-up was 26 months (range 22-30 months). The authors conclude that Tutopatch is a safe wrapping material.

- KROLL, J ET AL.  
 Defect covering of a perforated corneal ulcer with bovine pericardium transplant (Defektdeckung eines perforierten Hornhautulkus mit einem bovinen Perikardtransplantat).  
 Der Ophthalmologe, 111(1); 2014  
 An operative treatment of an 89-year-old patient presented in an emergency situation with a perforated corneal ulcer of neurotrophic genesis was performed. Tutopatch® represents a good alternative for covering perforated corneal ulcers. (Article in German)
- ORTUÑO-PRADOS VJ ET AL.  
 Treatment of a neurotrophic corneal ulcer with solid platelet-rich plasma and Tutopatch®  
 Arch Soc Esp Oftalmol, 86(4); 2011  
 Case report of a patient with neurotrophic queratitis in the left eye treated with a Tutopatch® cover and platelet-rich plasma (PRP). The authors found this form of treatment very effective for progressive ulcers. Tutopatch® may constitute an alternative to amniotic membrane transplantation.
- RINNA C ET AL.  
 Orbital floor restoration  
 J Craniofac Surg, 16(6); 2005  
 Retrospective study with 379 patients treated over a period of 9 years. Smaller fractures (<1cm diameter) were treated with a collagen membrane alone, larger fractures received various rigid materials covered with a collagen membrane. A total of 268 patients received a Tutopatch. No infection occurred and all cases of enophthalmus were corrected.
- SCHMAEL F ET AL.  
 Preseptal transconjunctival approach for orbital floor fracture repair: ophthalmologic results in 209 patients  
 Acta Oto-Laryngologica, 126(4); 2006  
 Retrospective study on the results of 209 orbital floor reconstructions via a special approach. Defects were covered with Tutopatch.
- SCUOLA A ET AL.  
 Bovine Pericardium for Scleral Closure in Transscleral Local Resection of Choroidal Melanoma  
 Retina, 28(10); 2008  
 Results of successful using bovine pericardium patches are presented in order to strengthen globe stability, repairing possible iatrogenic defects in the sclera flap, and to assist immediate vitrectomy, thus avoiding fluid leakage at the margins of the sutured scleral flap.
- YALNIZ-AKKAYA Z ET AL.  
 Late-Onset Glaucoma-Filtrating Bleb Leak in a Penetrating Keratoplasty Patient: A Case Report  
 Case Reports in Ophthalmological Medicine in: Hindawi; 2012  
 Surgical bleb revision using Tutopatch and conjunctival advancement in a 48-year-old female patient was performed. The authors conclude that this method seems to be effective for treating late bleb leaks. However, careful follow-up is required for detecting recurrent leaks and elevated intraocular pressure.  
<http://www.hindawi.com/journals/criopm/2012/810751/>
- YALVAC IS ET AL.  
 Double-layer pericardium sandwich technique of Ahmed glaucoma valve implantation in patient with anterior necrotizing scleritis  
 Techniques in Ophthalmology, 3(2); 2005  
 Case report. Tutopatch was used to reinforce the thinned sclera and for coverage of a glaucoma shunt.

- YOO C ET AL.  
Pericardium Plug in the Repair of the Corneoscleral Fistula After Ahmed Glaucoma Valve Explantation  
Korean Journal of Ophthalmology, (22); 2008  
A Tutopatch plug was used to repair a corneoscleral fistula after Ahmed Glaucoma Valve explantation in 4 cases. No complication related to the fistula developed at the latest follow-up (range: 12-26 months).

### **Ear Nose Throat**

- DE DORLODOT, C ET AL.  
Are bovine pericardium underlay xenograft and butterfly inlay autograft efficient for transcanal tympanoplasty?  
Eur Arch Otorhinolaryngol, Dec; 2013  
The authors recommend the butterfly technique for non-marginal perforation of the tympanic membrane after excision of the perforation edge and Tutopatch for bigger perforation or when standard autografts are not available
- BAUMANN I ET AL.  
Ossifying fibroma of the ethmoid involving the orbit and the skull base  
Otolaryngology-Head and Neck Surgery, 133; 2005  
Case report. The skull base was reconstructed with Tutopatch.
- CAVAZZA S ET AL.  
Two anomalous localizations of mucocoele: clinical presentation and retrospective review  
Acta Otorhinolaryngologica Italica, 27; 2007  
Two case reports. In one a defect in the posterior wall of the frontal sinus was closed with Tutopatch.
- GERARD JM AND GERSDORFF M  
The Tutopatch graft for transcanal myringoplasty  
B-ENT, 2(4); 2006  
In 11 patients a transcanal myringoplasty was performed. Follow-up was 3-14 months. In 10 patients the ear drum could be successfully closed. Except one all patients had an air-bone gap of less than 10 dB at 0.5, 1, 2 and 4 kHz. The authors conclude that with Tutopatch the same results can be obtained as with autologous temporalis fascia, perichondrium or cartilage-perichondrium.
- INCESULU A ET AL.  
The use of autograft and xenograft material in tympanoplasty: short term results  
KBB-Forum, 1(1); 2002  
Results of 31 patients with autologous temporalis fascia and 29 with Tutopatch for tympanoplasty were compared. Graft survival rate was 78% in the autograft group and 73% in the Tutoplast group.  
[http://kbb-forum.net/journal/summary\\_en.php3?id=7](http://kbb-forum.net/journal/summary_en.php3?id=7)
- Michel O  
Transnasal surgery of the orbita (Transnasale Chirurgie der Orbita)  
HNO, 48; 2000  
Overview. For reconstruction of the medial orbital wall the author recommends Tutopatch, Tutoplast Fascia lata or Tutoplast Fascia temporalis. He advises against a polydioxanon foil because of difficult handling due to its rigidity and very slow epithelisation. (Article in German)

- RINNA C ET AL.  
Orbital floor restoration  
J Craniofac Surg, 16(6); 2005  
Retrospective study on 379 patients. 268 had Tutopatch. The authors discuss the pros and cons of synthetic materials.
- SCHMAEL F ET AL.  
Preseptal transconjunctival approach for orbital floor fracture repair: ophthalmologic results in 209 patients  
Acta Oto-Laryngologica, 126; 2006  
Retrospective study. 48 patients received Tutopatch to cover the defect. No complications attributable to the graft were recorded. 1 PDS patch dislocated.
- SCHNEIDER G  
Bioimplants-characteristics and use (Bioimplantate-Eigenschaften und Anwendungshinweise)  
Laryngo Rhino Otol, 83; 2003  
CME article describing various biological grafts in ENT. Tutopatch is recommended for duraplasty, tympanoplasty and septoplasty. (Article in German, English translation available)

### **Cardiothoracic**

- US MH ET AL.  
A retrospective comparison of bovine pericardium and polytetrapolyethylene patch for closure of ventricular septal defects  
J Int Med Res, 32; 2004  
The outcome of 22 patients treated with Tutopatch and 57 treated with PTFE were compared. The outcome in the Tutopatch group was uneventful in all patients while in the PTFE group a patch dehiscence occurred in 4 patients. The authors prefer Tutopatch because of its handling characteristics, elasticity and the lower risk of endocarditis.
- US MH ET AL.  
Surgical correction in a patient with homozygous familial hypercholesterolemia  
Turkish J Thorac Cardiovasc Surg, 13(4); 2005  
Case report. The stenotic aortic root was enlarged with Tutopatch.
- BASLAIM G AND BASIONI A  
Repair of complete atrioventricular septal defects. Results with maintenance of the coronary sinus on the right atrial side  
J Card Surg, 21; 2006  
Retrospective study on 51 patients. Autologous pericardium was used for 25 ventricular septal defects and 40 atrial septal defects. Bovine pericardium (Tutopatch or Supple Peri-Guard) was used in 13 ventricular and 10 atrial septal defects.

### **Abdominal**

- SIMON S ET AL.  
Experience with prosthetic implants in congenital abdominal wall defects  
Abstract. German Association of Surgeons, annual Meeting, 2001  
Overview over 65 patients treated between 1984 and 2000. Until 1998 Tutoplast Dura was used and Tutopatch thereafter. The authors look at Tutopatch as an appropriate material because of the absence of immunological rejections and total graft transformation into connective tissue.

- SAXENA AK ET AL.  
Delayed three-stage closure of giant omphalocele using pericard patch  
Hernia, 12; 2008  
Case report of a giant omphalocele that required a three-stage delayed closure and was managed using Tutopatch. Delayed staged closure has the advantage of avoiding multiple surgical procedures in the neonatal period and reduces morbidity in the management of giant omphaloceles.
- VAN TUIL C ET AL.  
Experience with management of anterior abdominal wall defects using bovine pericard  
Hernia, 10; 2006  
Report on 24 patients with congenital abdominal wall defects treated with Tutopatch over a period of 5 years. The authors conclude that Tutopatch is an optimal biomaterial for the closure of anterior abdominal wall defects in gastrochisis and omphalocele
- YAVUZ N ET AL.  
Laparoscopic repair of Morgagni hernia  
Surg Laparosc Endosc Percutan Tech, 16(3); 2006  
Report on the successful treatment of a Morgagni hernia in 5 patients. In one patient Tutopatch was used.

### **Breast Reconstruction**

- SEMPRINI G ET AL.  
The bovine pericardial patch in breast reconstruction: a case report.  
Il Giornale di chirurgia, 33(11/12); 2012  
Tutopatch® has been used in immediate breast reconstruction after mastectomy to close laterally the subpectoral pocket, allowing bigger volume prosthesis to be placed. There were no particular postoperative complications, and after 12 months of follow up the authors have found a valid functional and aesthetic result. The authors consider Tutopatch® as a valid alternative to other acellular dermal matrices specifically designed for breast reconstruction.  
<http://www.giornalechirurgia.it/common/php/portiere.php?ID=88f319b1cb6b1580fa19f1573d707fe76>

### **Animal studies**

- ASLAN M ET AL.  
Guided bone regeneration (GBR) on healing bone defects: A histological study in rabbits  
J Contemp Dent Practice, 5(2); 2004  
Two 3mm bone defects were drilled in the right tibia of 30 rabbits. One was covered with Tutopatch and the other left empty as control. 10 rabbits each were sacrificed on day 10, 20 and 30 and the defects examined histologically. In all Tutopatch covered defects new bone formation was found while in 4 control cases fibrous tissue only was seen. No resorption or dehiscence of Tutopatch occurred.
- AYYILDIZ A ET AL.  
A comparison of free skin graft, fascia lata, alloderm, bovine pericardium and primary repair in urethrocutaneous fistulas without diversion: an experimental study  
Pediatr Surg Int, 22; 2006  
Comparison of several collagen membranes, among them Tutopatch (not specifically named but confirmed by the corresponding author via e-mail). Various parameters were rated with points. Tutopatch reached the highest number.

- ÇÖLOĞLU H ET AL.  
Comparison of autogenous cartilage, acellular dermis, and solvent-dehydrated pericardium for the prevention and correction of dorsal nasal irregularities: an experimental study  
Aesth Plast Surg, 36; 2012  
30 New Zealand rabbits underwent rhinoplasty with different materials among them with Tutopatch. The major advantages of Tutopatch were the ease of obtaining it without any donor site morbidity, shorter operative procedures, and lower distortion rates due to lack of cartilage memory. There were no adverse effects such as ulceration or extrusion and no adhesion to the overlying skin. In contrast, the used autogenous ear cartilage caused more cartilage resorption than expected by the authors.  
The authors conclude that Tutopatch may be used successfully as an "onlay" graft for dorsal nasal problems compared to autogenous cartilage, which is commonly used for this purpose.
- DULAURENT T ET AL.  
Use of bovine pericardium (Tutopatch®) graft for surgical repair of deep melting corneal ulcers in dogs and corneal sequestra in cats  
Veterinary Ophthalmology, 17(2); 2014  
The efficacy of Tutopatch in the treatment of deep melting corneal ulcers was evaluated in three dogs and corneal sequestra in three cats. Beginning corneal neovascularization was observed in all cases 1 week postoperatively. After 4 weeks, in two dogs and in all cats, the vascularization was regressing and the graft was integrated into the cornea, which was regaining transparency. Two months after the surgery, 5 of 6 corneas in two dogs and three cats had healed with focal corneal scarring. The authors conclude that Tutopatch offers a promising option for surgical reconstruction of the cornea following keratectomy for the management of corneal ulcers and sequestra.
- KAPAN S ET AL.  
Comparison of PTFE, pericardium bovine and fascia lata for repair of incisional hernia in rat model, experimental study  
Hernia, 7(1); 2003  
A 1x2cm subcutaneous full thickness abdominal wall defect was created in rats and closed with patches of Tutopatch, Tutoplast Fascia lata and Goretex Dual Mesh. After one and two weeks bursting pressure, tensile strength of the graft-tissue connection, histological graft maturation and adhesion formation were evaluated. All values were comparable for Tutopatch and GoreTex while bursting strength was lower with Fascia lata. 7
- KOHL T ET AL  
Percutaneous fetoscopic patch coverage of experimental lumbosacral full-thickness skin lesions in sheep – a minimally invasive technique to minimize maternal trauma from fetal surgery for myelomeningocele  
Surg Endosc, 17(8); 2003  
In 7 sheep fetuses a full thickness defect in the lumbosacral region was created endoscopically to mimic a fetal myelomeningocele. In 6 the defect was closed with e-PTFE and in one with Tutopatch. After delivery 45 to 55 days later an autopsy was performed. The e-PTFE patches were still present and could be easily removed while the Tutopatch was completely remodelled.
- ROTHAMEL D ET AL.  
Biodegradation of differently crosslinked collagen membranes: an experimental study in the rat  
Clin Oral Implants Res, 16; 2005  
Various collagen membranes, among them Tutopatch, were implanted into dorsal muscle pouches in rats and examined histologically and histomorphometrically after 2, 4, 8, 16 24 weeks. Tutopatch was completely vascularised after 4 weeks and nearly completely remodelled after 16 weeks.

- SCHWARZ F ET AL.  
 Angiogenesis pattern of native and cross-linked collagen membranes: an immunohistochemical study in the rat  
 Clin Oral Implants Res, 17(4); 2006  
 Various commercially available collagen membranes, among them the Tutodent membrane, and 3 experimental collagen membranes were tested. 4 each, randomly selected, were implanted into muscle pouches in the back of rats and examined after 2, 4, 8, 16 and 24 weeks. The Tutodent membrane showed a homogenous transmembraneous formation of blood vessels 4-8 weeks after implantation.
- SUNGUR N ET AL.  
 Prevention of tendon adhesions by the reconstruction of the tendon sheet with solvent dehydrated bovine pericard. An experimental study.  
 J Trauma, 61(6); 2006  
 Chicken were chosen as a model because of the anatomical and functional similarity to the human flexor tendon apparatus. In 30 chicken the tendon sheet between the proximal and distal ring band was resected. In 10 animals the defect remained untreated, in 10 it was covered with autologous fascia and in 10 with Tutopatch. Animals were sacrificed after 5 and 12 weeks. The gliding function in the Tutopatch group was significantly better than in the two other groups.
- QIU QQ ET AL.  
 In Vitro and In Vivo Comparison of Tutopatch® Bovine Pericardium and Veritas® Collagen Matrix for Hernia Repair  
 White Paper, RTI Surgical, Inc., Alachua; FL, USA; 2014  
 Nine dogs received Tutopatch® on one side and Veritas on the contralateral side sutured over the abdominal wall. Tutopatch® showed more favourable characteristics at all investigation times (4, 8 and 12 weeks) postop. As Example Tutopatch® did not contract and remained tension free in vivo and exhibited significantly greater cell infiltration and neovascularization than Veritas throughout the 12-week implantation period.

### **Laboratory studies**

- ROTHAMEL D ET AL.  
 Biocompatibility of various collagen membranes in cultures of human PDL fibroblasts and human osteoblast-like cells  
 Clin Oral Implants Res, 15(4); 2004  
 Various collagen membranes, among them Tutopatch, were tested in cell cultures for cell proliferation. While cell count was highest in controls, cell shape was the same on Tutopatch than in controls. It was concluded that Tutopatch promotes cell attachment and proliferation.
- MODER D ET AL.  
 Collagenase digestion of collagen membranes in vitro  
 Abstract of a presentation, IADR-Continental European and Israeli Divisions, annual Meeting, 2007  
 Degradation kinetics of BioGuide, BioMend, Ossix, Tutodent and an experimental membrane after incubation in collagenase were analyzed. BioGuide was completely digested after 2h and BioMend and Ossix after 6h. After 12h 40% of Tutodent remained.  
[http://iadr.confex.com/iadr/israel07/techprogram/abstract\\_95362.htm](http://iadr.confex.com/iadr/israel07/techprogram/abstract_95362.htm)